



**Office of Controller General of Patents
Designs, and Trademarks**
Indian Patent Information Retrieval System



(IPIRS)

01 March 2010 17:28	Guidelines For Searching	Home
---------------------	--------------------------	------

**Published Patent
Applications**

(12) PATENT APPLICATION PUBLICATION

(21) Application No. : 701/DEL/2004

(19) INDIA

(22) Date of filing of Application :12/04/2004

(43) Publication Date : 16/06/2006

Journal No. - 24/2006

(54) Title of the invention : GENITALIA EXTERNAL PROTHESIS FOR TREATMENT OF MALE EXTERNAL GENITALIA HYPOPLASIA.

(51) International classification	:A61H 1/00
(31) Priority Document No	:EP0338101.0
(32) Priority Date	:28/04/2003
(33) Name of priority country	:Spain
(86) International Application No	:NA
Filing Date	:NA
(87) International Publication No	: NA
(61) Patent of Addition to	:NA
Application Number	:NA
Filing Date	
(62) Divisional to Application	:NA
Number	:NA
Filing Date	

(71)Name of Applicant :
1)GOMEZ DE DIEGO, EDUARDO
Address of Applicant :Gran Via 6, E-28013-Madrid(SE) Spain
(72)Name of Inventor :
1)GOMEZ DE DIEGO, EDUARDOGO (Spain)

(57) Abstract :

The invention describes a device designed for the treatment of the external male genitalia hipoplasia, creating a growth of tissues and a correction of existing curvatures. The invention's device is unique for using an external genitalia fixation(1;13) associated to a dynamometric system (5) and to a gland protection system(15). The system of external genitalia inferior fixation(1) and superior (13) produces on the male genitalia a traction force. The gland protection system (5) creates on the male genitalia a tissue expansion. The gland protection system(15) avoids possible cuts, this allows a long term use of the device, free of discomfort. The device's user can wear it for many months, during a few hours a day. This way, the chance of a tissue expansion plus an elongated expansion in time are fundamental for the external male genitalia hipoplasia treatment.

Number of Pages = 8